

E 494
Insert
for the gentle holding of 5 microtiter plates.



EAN: 4002516853947 / Material number: 12669920 /
Old Material Number: 69549402D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 6011	•
PLW 6111	•
PLW 7111	•
PLW 8505	•
Product category	
Other inserts	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of different utensils	•
Capacity	
Microtiter plates [number]	5
Dimensions and weight	
External dimensions, gross height in mm	90
External dimensions, gross width in mm	230
External dimensions, gross depth in mm	465
Net weight in kg	0.48
Gross weight in kg	0.72