


## A 613 Module

for the optimum reprocessing of volumetric and measuring pipettes.



- Equipped with 98 injector nozzles
- Suitable for measuring pipettes with a volume of 1–50 ml
- Suitable for volumetric pipettes with a volume of 1–100 ml
- For use in A 503 mobile unit

EAN: 4002516854128 / Material number: 12670020 / Old Material Number: 69661301D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 8615	•
PLW 8616	•
PLW 8617	•
<b>Product category</b>	
Insert/module for pipettes	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Area of application</b>	
Reprocessing of pipettes/phials	•
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Pipettes [number]	98
<b>Dimensions and weight</b>	
External dimensions, net height in mm	194
External dimensions, net width in mm	254
External dimensions, net depth in mm	604
External dimensions, gross height in mm	217
External dimensions, gross width in mm	283
External dimensions, gross depth in mm	615
Net weight in kg	2.69
Gross weight in kg	3.05