

E 496

Injector nozzle

with retaining clamp for intubation material.



EAN: 4002516854579 / Material number: 12670260 / Old Material Number: 69749602D

Product category

Washer-disinfectors, medical

•

Large-capacity washer-disinfectors, medical

•

Machine affiliation

G 7823

•

G 7824

•

G 7825

•

G 7826

•

PG 8527

•

PG 8528

•

PG 8535

•

PG 8536

•

PG 8562

•

Product category

Injector nozzle

•

Product properties

Material

Stainless steel, Plastic

Colour

Stainless Steel, Black

Area of application

Reprocessing of anaesthetic instruments

•

Reprocessing of MIS instruments

•

Reprocessing of lumened instruments

•

Dimensions and weight

External dimensions, net height in mm

120

External dimensions, gross height in mm

12

External dimensions, gross width in mm

80

External dimensions, gross depth in mm

120

Net weight in kg

0.01

Gross weight in kg

0.01