

E 792

Silicone hose with connection
for the connection to female Luer lock connectors.



- With male Luer lock connection
- Length 160 mm
- 4 silicone hoses per pack

EAN: 4002516854760 / Material number: 12670280 / Old Material Number: 69779202D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Other accessories	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Transparent
Area of application	
Reprocessing of MIS instruments	•
Reprocessing of ophthalmology instruments	•
Reprocessing of lumened instruments	•
Dimensions and weight	
External dimensions, gross height in mm	12
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	120
Net weight in kg	0.08
Gross weight in kg	0.09