


E 803

Insert

for the optimum loading of up to 30 ear and nose specula.



- 160 compartments approx. 13 x 13 mm
- Mesh width on base: 1.7 mm
- For use in upper and lower basket

EAN: 4002516854715 / Material number: 12670290 / Old Material Number: 69580302D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8534	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Insert for ENT instruments	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	160
Area of application	
Reprocessing of ENT instruments	•
Dimensions and weight	
External dimensions, gross height in mm	132
External dimensions, gross width in mm	233
External dimensions, gross depth in mm	484
Net weight in kg	0.80
Gross weight in kg	1.07