

A 621

Injector module

for the optimum reprocessing of laboratory glassware, 20 injector positions.



- Can be equipped with up to 20 nozzles based on individual needs
- For use in A 503 and APLW 511 mobile units
- Part of the EasyLoad system

EAN: 4002516855125 / Material number: 12678110 / Old Material Number: 69662101D

Product category	
Large-chamber laboratory washers	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Insert/module for narrow-necked glassware	•
Insert for various flasks	•
Area of application	
Reprocessing of laboratory glassware	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Part of the EasyLoad system	•
Capacity	
Narrow-necked glassware [number]	20
Dimensions and weight	
External dimensions, net height in mm	113
External dimensions, net width in mm	173
External dimensions, net depth in mm	608
External dimensions, gross height in mm	170
External dimensions, gross width in mm	208
External dimensions, gross depth in mm	645
Net weight in kg	1.93
Gross weight in kg	2.27