

A 621

Injector module

for the optimum reprocessing of laboratory glassware, 20 injector positions.



- Can be equipped with up to 20 nozzles based on individual needs
- For use in A 503 mobile unit
- Part of the EasyLoad system

EAN: 4002516855125 / Material number: 12678110 / Old Material Number: 69662101D

| | |
|--|-----------------|
| Product category | |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| PLW 8615 | • |
| PLW 8616 | • |
| PLW 8617 | • |
| Product category | |
| Insert/module for narrow-necked glassware | • |
| Insert for various flasks | • |
| Product properties | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Part of the EasyLoad system | • |
| Area of application | |
| Reprocessing of laboratory glassware | • |
| Capacity | |
| Narrow-necked glassware [number] | 20 |
| Dimensions and weight | |
| External dimensions, net height in mm | 113 |
| External dimensions, net width in mm | 173 |
| External dimensions, net depth in mm | 608 |
| External dimensions, gross height in mm | 170 |
| External dimensions, gross width in mm | 208 |
| External dimensions, gross depth in mm | 645 |
| Net weight in kg | 1.93 |
| Gross weight in kg | 2.27 |