

E 936

Insert

for optimum loading of up to 15 lamp handles.



EAN: 4002516854326 / Material number: 12678190 / Old Material Number: 69593602D

## Product category

Large-capacity washer-disinfectors, medical

## Machine affiliation

G 7823

G 7824

G 7826

PG 8527

PG 8528

## Product category

Other inserts

## Area of application

Reprocessing of surgical instruments

## Product properties

Material

Stainless steel

Colour

Stainless Steel

## Capacity

Lamp handles [number]

15

## Dimensions and weight

External dimensions, net height in mm

90

External dimensions, net width in mm

245

External dimensions, net depth in mm

500

External dimensions, gross height in mm

135

External dimensions, gross width in mm

277

External dimensions, gross depth in mm

529

Net weight in kg

0.64

Gross weight in kg

1.05