

E 975/2

Mobile unit

for the optimum cleaning of instruments on 2 levels.



- Equipped with an integrated spray arm
- Magnetic strip for automatic mobile unit recognition
- The vertical clearance for both levels is different

EAN: 4002516853916 / Material number: 12678470 / Old Material Number: 69597504D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Large-chamber washer-disinfectors, fire brigades and emergency services	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Mobile unit for inserts	•
Mobile unit for mesh trays	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of surgical instruments	•
Capacity	
DIN mesh trays [number]	6
SPRI mesh trays [number]	4
Dimensions and weight	
External dimensions, net height in mm	426
External dimensions, net width in mm	636
External dimensions, net depth in mm	790
External dimensions, gross height in mm	790
External dimensions, gross width in mm	715
External dimensions, gross depth in mm	860
Net weight in kg	11.40
Gross weight in kg	28.60
Companion pack	
Magnetic strip for automatic mobile unit recognition (incl. magnets) [number]	1