

A 602

Module

For the optimum reprocessing of instruments used in robot-assisted surgery.



- Suitable for 8 instruments, latest generation (Xi)

EAN: 4002516843108 / Material number: 12678550 / Old Material Number: 69660202D

Product category

Large-capacity washer-disinfectors, medical

•

Machine affiliation

PG 8527

•

PG 8528

•

Product category

Module/rack for robotic instruments

•

Area of application

Reprocessing of robotic instruments

•

Product properties

Material

Stainless steel

Colour

Stainless Steel

Capacity

Robotics instruments [number]

8

Dimensions and weight

External dimensions, net height in mm

303

External dimensions, net width in mm

555

External dimensions, net depth in mm

318

External dimensions, gross height in mm

353

External dimensions, gross width in mm

354

External dimensions, gross depth in mm

590

Net weight in kg

5.15

Gross weight in kg

6.39