

A 805 Holder

For the optimum reprocessing of instruments used in robot-assisted surgery.



- For special instruments, surgical clamps (Si)

EAN: 4002516854272 / Material number: 12680080 /
Old Material Number: 69680502D

Product category	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Module/rack for robotic instruments	•
Area of application	
Reprocessing of robotic instruments	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Capacity	
Robotics instruments [number]	1
Dimensions and weight	
External dimensions, net height in mm	92
External dimensions, net width in mm	41
External dimensions, net depth in mm	165
External dimensions, gross height in mm	73
External dimensions, gross width in mm	163
External dimensions, gross depth in mm	245
Net weight in kg	0.19
Gross weight in kg	0.30