

A 806

Holder

For the optimum reprocessing of instruments used in robot-assisted surgery.



- For special instruments with side rinse connection

EAN: 4002516854289 / Material number: 12680100 /  
Old Material Number: 69680602D

<b>Product category</b>	
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8527	•
PG 8528	•
<b>Product category</b>	
Module/rack for robotic instruments	•
<b>Area of application</b>	
Reprocessing of robotic instruments	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Capacity</b>	
Robotics instruments [number]	1
<b>Dimensions and weight</b>	
External dimensions, net height in mm	77
External dimensions, net width in mm	41
External dimensions, net depth in mm	165
External dimensions, gross height in mm	70
External dimensions, gross width in mm	160
External dimensions, gross depth in mm	240
Net weight in kg	0.16
Gross weight in kg	0.27