


## A 807 Holder

For the optimum reprocessing of instruments used in robot-assisted surgery.

EAN: 4002516854685 / Material number: 12680120 / Old Material Number: 69680702D



- For special instruments, surgical clamps, older generation, S/Si

<b>Product category</b>	
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8527	•
PG 8528	•
<b>Product category</b>	
Module/rack for robotic instruments	•
<b>Area of application</b>	
Reprocessing of robotic instruments	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Capacity</b>	
Robotics instruments [number]	1
<b>Dimensions and weight</b>	
External dimensions, net height in mm	77
External dimensions, net width in mm	41
External dimensions, net depth in mm	165
External dimensions, gross height in mm	25
External dimensions, gross width in mm	245
External dimensions, gross depth in mm	485
Net weight in kg	0.16
Gross weight in kg	0.24