

## A 808 Holder

For the optimum reprocessing of instruments used in robot-assisted surgery.



- For special instruments, surgical clamps (Xi)

EAN: 4002516853886 / Material number: 12680130 /  
Old Material Number: 69680802D

### Product category

Large-capacity washer-disinfectors, medical

•

### Machine affiliation

PG 8527

•

PG 8528

•

### Product category

Module/rack for robotic instruments

•

### Area of application

Reprocessing of robotic instruments

•

### Product properties

Material

Stainless steel

Colour

Stainless Steel

### Capacity

Robotics instruments [number]

1

### Dimensions and weight

External dimensions, net height in mm

50

External dimensions, net width in mm

145

External dimensions, net depth in mm

159

External dimensions, gross height in mm

76

External dimensions, gross width in mm

164

External dimensions, gross depth in mm

244

Net weight in kg

0.21

Gross weight in kg

0.31