

MF 27/28-1

Transfer trolley

for loading & unloading washer-&disinfectors.



- Equipped with 4 castors, 2 of which can be locked
- Simple and convenient handling thanks to rounded handles
- Dockable or lockable to rolling conveyor etc. with both sides
- Infeed height 850 mm, -100 mm, +150 mm
- With removable drip tray H 70, W 603, D 866 mm

EAN: 4002516853923 / Material number: 12680160 / Old Material Number: 69200112D

Product category

Large-capacity washer-disinfectors, laboratory

•

Large-capacity washer-disinfectors, medical

•

Machine affiliation

PG 8527

•

PG 8528

•

Product category

Storage/transport system

•

Product properties

Material

Stainless steel

Colour

Stainless Steel

Dimensions and weight

External dimensions, net height in mm

1050

External dimensions, net width in mm

740

External dimensions, net depth in mm

930

External dimensions, gross height in mm

1155

External dimensions, gross width in mm

800

External dimensions, gross depth in mm

1200

Net weight in kg

28.00

Gross weight in kg

47.60