

MF/3

Trolley

for loading & unloading washer-&disinfectors.



- For easy handling of the mobile instrument unit
- With lifting technology operated by foot pressure
- Equipped with 2 lockable castors
- Infeed height 751 mm, +/- 100 mm

EAN: 4002516853817 / Material number: 12680170 / Old Material Number: 69200119D

## Product category

Washer-disinfectors, laboratory

•

Washer-disinfectors, medical

•

## Machine affiliation

G 7823

•

G 7824

•

G 7825

•

G 7826

•

## Product category

Storage/transport system

•

## Product properties

Material

Stainless steel

Colour

Stainless Steel

## Dimensions and weight

External dimensions, net height in mm

1182

External dimensions, net width in mm

660

External dimensions, net depth in mm

807

External dimensions, gross height in mm

770

External dimensions, gross width in mm

800

External dimensions, gross depth in mm

1200

Net weight in kg

34.20

Gross weight in kg

55.80