

UGT 44-60/85

Box plinth, with door

for ergonomic loading and unloading of the machine.



- Offers the option of storing e.g. process chemicals

EAN: 4002516854906 / Material number: 12680260 /
Old Material Number: 69506002D

Product category

Washer-disinfectors, laboratory

•

Thermal disinfectors

•

Washer-disinfectors, medical

•

Machine affiliation

PG 8581

•

PG 8582

•

PG 8583

•

PG 8591

•

PG 8592

•

PG 8593

•

Product category

Plinth/base

•

Product properties

Material

Stainless steel

Colour

Stainless Steel

Dimensions and weight

External dimensions, gross height in mm

468

External dimensions, gross width in mm

568

External dimensions, gross depth in mm

610

Net weight in kg

18.20

Gross weight in kg

19.80