


RP 900

Plinth, on wheels

For machines 900 mm wide, incl. drip tray without openings for floor drain.



- Suitable for single-door machines, 900 mm wide
- Castors allow the machine to be pulled forward for servicing
- Only suitable for machines with a drain pump
- Easy access for servicing

EAN: 4002516855170 / Material number: 12680650 / Old Material Number: 69253012D

Product category

Large-capacity washer-disinfectors, laboratory

Machine affiliation

PLW 8615

Product category

Plinth/base

Product properties

Material

Stainless steel

Colour

Stainless Steel

Dimensions and weight

External dimensions, net height in mm

100

External dimensions, net width in mm

900

External dimensions, net depth in mm

990

External dimensions, gross height in mm

155

External dimensions, gross width in mm

940

External dimensions, gross depth in mm

1200

Net weight in kg

30.80

Gross weight in kg

33.10