

TC 1150

Dresser panelling

For machines 1150 mm wide without a steam condenser or drying unit.



- Includes service door with lock
- Diaphragm wall panelling up to ceiling must be provided on site

EAN: 4002516853848 / Material number: 12680670 /
Old Material Number: 69252005D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8617	•
Product category	
Cover	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Dimensions and weight	
External dimensions, net height in mm	1002
External dimensions, net width in mm	835
External dimensions, net depth in mm	1200
External dimensions, gross height in mm	1002
External dimensions, gross width in mm	860
External dimensions, gross depth in mm	1200
Net weight in kg	60.13
Gross weight in kg	75.13