

VE P 2000/1
Full demineralisation cartridge
For optimum water preparation.



- Max. delivery: 450 l/h
- ConductivityMeter CM/1 measures remaining capacity
- Max. water inlet temperature 40 °C – pressure range from 0.5–10 bar
- Height 40 cm, Ø 24 cm

EAN: 4002516874010 / Material number: 12680820 / Old Material Number: 69431102D

Product category	
Fresh-water dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
G 8066	•
PG 8067	•
PG 8041	•
PG 8043	•
PG 8045	•
PG 8055	•
PG 8055 U	•
PG 8056	•
PG 8056 U	•
PG 8057 TD	•
PG 8057 TD U	•
PG 8058	•
PG 8058 U	•
PG 8059	•
PG 8059 U	•
PG 8060	•
PG 8061	•
PG 8061 U	•
PG 8063	•
PG 8164	•
PG 8165	•
PG 8166	•
PG 8169	•
PG 8172	•
PG 8172 Eco	•
PTD 901	•
PG 8096	•
PG 8096 U	•
PG 8099	•
PG 8099 U	•
PFD 400	•
PFD 400 U	•
PFD 401	•
PFD 401 U	•
PFD 402	•
PFD 402 U	•
PFD 404	•
PFD 404 U	•
PFD 405	•
PFD 405 U	•

VE P 2000/1
Full demineralisation cartridge
For optimum water preparation.



EAN: 4002516874010 / Material number: 12680820 / Old Material Number: 69431102D

PFD 407	•
PFD 407 U	•
PFD 408 U	•
PTD 701	•
PTD 702	•
PTD 703	•
PTD 704	•
PG 8581	•
PG 8582	•
PG 8582 CD	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
Product category	
Water treatment	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Dimensions and weight	
External dimensions, gross height in mm	250
External dimensions, gross width in mm	260
External dimensions, gross depth in mm	610
Net weight in kg	11.87
Gross weight in kg	12.60