

MC/2

Transfer trolley

height adjustable, for loading & unloading washer-&disinfectors.



- Equipped with 4 castors, 2 of which can be locked
- 2 loading racks (trough shaped) as well as handle and cover plate

EAN: 4002516853855 / Material number: 12680970 / Old Material Number: 69200113D

## Product category

Washer-disinfectors, medical

•

## Machine affiliation

PG 8535

•

PG 8536

•

PG 8581

•

PG 8582

•

PG 8582 CD

•

PG 8583

•

PG 8583 CD

•

PG 8591

•

PG 8592

•

PG 8593

•

## Product category

Storage/transport system

•

## Product properties

Material

Stainless steel

Colour

Stainless Steel

## Dimensions and weight

External dimensions, net height in mm

1003

External dimensions, net width in mm

600

External dimensions, net depth in mm

805

External dimensions, gross height in mm

1210

External dimensions, gross width in mm

800

External dimensions, gross depth in mm

1200

Net weight in kg

18.00

Gross weight in kg

47.80