

MT  
Trolleys  
for loading & unloading washer-&disinfectors.



- For the storage and transport of baskets and inserts
- Divided into 4 levels that can be adjusted by 102.5 mm
- Moveable with 4 lockable castors
- Loading volume 549 x 599 mm

EAN: 4002516854159 / Material number: 12681000 / Old Material Number: 69200114D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8582	•
PG 8582 CD	•
PG 8583	•
PG 8583 CD	•
PG 8592	•
PG 8593	•
<b>Product category</b>	
Storage/transport system	•
<b>Dimensions and weight</b>	
External dimensions, gross height in mm	920
External dimensions, gross width in mm	775
External dimensions, gross depth in mm	2001
Net weight in kg	23.00
Gross weight in kg	37.60