

TE P 2800/1
 Partial demineralisation cartridge
 For optimum water preparation.



- Max. delivery: 800 l/h
- FlowMeter FM measures remaining capacity
- Max. water inlet temperature 40 °C – pressure range from 0.5–10 bar
- Height 56 cm, Ø 24 cm
- Softening capacity at carbonate hardness of 10 °dH: 18 m³

EAN: 4002516873945 / Material number: 12681030 / Old Material Number: 67809502D

Product category	
Fresh-water dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Machine affiliation	
G 8066	•
PG 8067	•
PG 8041	•
PG 8043	•
PG 8045	•
PG 8055	•
PG 8055 U	•
PG 8056	•
PG 8056 U	•
PG 8057 TD	•
PG 8057 TD U	•
PG 8058	•
PG 8058 U	•
PG 8059	•
PG 8059 U	•
PG 8060	•
PG 8061	•
PG 8061 U	•
PG 8063	•
PG 8164	•
PG 8165	•
PG 8166	•
PG 8169	•
PG 8172	•
PG 8172 Eco	•
PTD 901	•
PG 8096	•
PG 8096 U	•
PG 8099	•
PG 8099 U	•
PFD 400	•
PFD 400 U	•
PFD 401	•
PFD 401 U	•
PFD 402	•
PFD 402 U	•
PFD 404	•
PFD 404 U	•
PFD 405	•
PFD 405 U	•
PFD 407	•
PFD 407 U	•
PFD 408 U	•

TE P 2800/1
Partial demineralisation cartridge
For optimum water preparation.

EAN: 4002516873945 / Material number: 12681030 / Old Material Number: 67809502D



PTD 701	•
PTD 702	•
PTD 703	•
PTD 704	•
Product category	
Water treatment	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black