

TE P 2800/1
 Partial demineralisation cartridge
 For optimum water preparation.



- -
- -
- -
- -
- -

EAN: 4002516873945 / Material number: 12681030 / Old Material Number: 67809502D

Product category	
Fresh-water dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Machine affiliation	
G 8066	•
PG 8067	•
PG 8041	•
PG 8043	•
PG 8045	•
PG 8055	•
PG 8055 U	•
PG 8056	•
PG 8056 U	•
PG 8057 TD	•
PG 8057 TD U	•
PG 8058	•
PG 8058 U	•
PG 8059	•
PG 8059 U	•
PG 8060	•
PG 8061	•
PG 8061 U	•
PG 8063	•
PG 8164	•
PG 8165	•
PG 8166	•
PG 8169	•
PG 8172	•
PG 8172 Eco	•
PTD 901	•
PG 8096	•
PG 8096 U	•
PG 8099	•
PG 8099 U	•
PFD 400	•
PFD 400 U	•
PFD 401	•
PFD 401 U	•
PFD 402	•
PFD 402 U	•
PFD 404	•
PFD 404 U	•
PFD 405	•
PFD 405 U	•
PFD 407	•
PFD 407 U	•
PFD 408 U	•

TE P 2800/1
Partial demineralisation cartridge
For optimum water preparation.



EAN: 4002516873945 / Material number: 12681030 / Old Material Number: 67809502D

PTD 701	•
PTD 702	•
PTD 703	•
PTD 704	•
Product category	
Water treatment	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Dimensions and weight	
External dimensions, net height in mm	240
External dimensions, net width in mm	235
External dimensions, net depth in mm	560
External dimensions, gross height in mm	250
External dimensions, gross width in mm	250
External dimensions, gross depth in mm	620
Net weight in kg	17.14
Gross weight in kg	17.80