

A 307

Insert

for the optimum loading of 20 theatre shoes up to size 48.



- Suitable for 20 theatre shoes
- Each with 265 mm tall supports
- For use in the lower basket carrier A 151
- Reprocessing of theatre shoes to size 48 possible

EAN: 4002516854975 / Material number: 12683150 / Old Material Number: 69630701D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
<b>Product category</b>	
Insert for theatre shoes	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of theatre shoes	•
<b>Capacity</b>	
Theatre shoes [number]	20
<b>Dimensions and weight</b>	
External dimensions, net height in mm	279
External dimensions, net width in mm	390
External dimensions, net depth in mm	360
External dimensions, gross height in mm	290
External dimensions, gross width in mm	396
External dimensions, gross depth in mm	480
Net weight in kg	1.52
Gross weight in kg	2.01