

A 308

Insert

For the optimum loading of 34 insoles for theatre shoes.



- Suitable for up to 34 insoles
- For use in upper or lower basket carriers

EAN: 4002516855415 / Material number: 12683170 / Old Material Number: 69630801D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
<b>Product category</b>	
Insert for theatre shoes	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of theatre shoes	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	71
External dimensions, net width in mm	475
External dimensions, net depth in mm	445
External dimensions, gross height in mm	270
External dimensions, gross width in mm	525
External dimensions, gross depth in mm	540
Net weight in kg	1.14
Gross weight in kg	2.00