

A 310
Insert
for the optimum loading of 8 theatre shoes.



- For use in upper baskets
- For loading 8 theatre shoes up to size 41
- Can be used in upper baskets A 101 and A 102

EAN: 4002516854661 / Material number: 12684420 / Old Material Number: 69631001D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Insert for theatre shoes	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of theatre shoes	•
Capacity	
Theatre shoes [number]	8
Dimensions and weight	
External dimensions, gross height in mm	100
External dimensions, gross width in mm	425
External dimensions, gross depth in mm	450
Net weight in kg	0.49
Gross weight in kg	0.83