

LP 2800/1

Demineralisation cartridge, empty

For the manual filling of single use resins.



- Max. delivery: 800 l/h
- Max. water inlet temperature 40 °C – pressure range from 0.5–10 bar
- Height 56 cm, Ø 24 cm

EAN: 4002516881124 / Material number: 12684720 / Old Material Number: 69431806D

Product category	
Fresh-water dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Laboratory washers	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
G 8066	•
PG 8067	•
PG 8041	•
PG 8043	•
PG 8045	•
PG 8055	•
PG 8055 U	•
PG 8056	•
PG 8056 U	•
PG 8057 TD	•
PG 8057 TD U	•
PG 8058	•
PG 8058 U	•
PG 8059	•
PG 8059 U	•
PG 8060	•
PG 8061	•
PG 8061 U	•
PG 8063	•
PG 8164	•
PG 8165	•
PG 8166	•
PG 8169	•
PG 8172	•
PG 8172 Eco	•
PTD 901	•
PG 8096	•
PG 8096 U	•
PG 8099	•
PG 8099 U	•
PFD 400	•
PFD 400 U	•
PFD 401	•
PFD 401 U	•
PFD 402	•
PFD 402 U	•
PFD 404	•
PFD 404 U	•
PFD 405	•
PFD 405 U	•

LP 2800/1

Demineralisation cartridge, empty

For the manual filling of single use resins.



EAN: 4002516881124 / Material number: 12684720 / Old Material Number: 69431806D

PFD 407	•
PFD 407 U	•
PFD 408 U	•
PTD 701	•
PTD 702	•
PTD 703	•
PTD 704	•
PG 8581	•
PG 8582	•
PG 8583	•
PG 8583 CD	•
PG 8591	•
PG 8592	•
PG 8593	•
Product category	
Water treatment	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel
Dimensions and weight	
External dimensions, net height in mm	240
External dimensions, net width in mm	240
External dimensions, net depth in mm	559
External dimensions, gross height in mm	620
External dimensions, gross width in mm	250
External dimensions, gross depth in mm	250
Net weight in kg	4.46
Gross weight in kg	5.12