

A 503

Mobile unit

no levels, for cleaning laboratory glassware/utensils.



- For loading modules and inserts
- Depending on the height, up to 6 modules can be used
- Upper rail and upper strut can be removed

EAN: 4002516854784 / Material number: 12685380 /
Old Material Number: 69650302D

Product category

Large-capacity washer-disinfectors, laboratory •

Machine affiliation

PLW 8615 •

PLW 8616 •

PLW 8617 •

Product category

Mobile unit with module connection •

Product properties

Material

Stainless steel, Plastic

Colour

Stainless Steel, Black

Area of application

Reprocessing of laboratory glassware •

Dimensions and weight

External dimensions, gross height in mm

790

External dimensions, gross width in mm

715

External dimensions, gross depth in mm

860

Net weight in kg

11.40

Gross weight in kg

22.20