

E 491
Insert
for the optimum loading of rotating instruments.



- Suitable for 18 turbine shafts
- Suitable for 30 contra angle shafts
- Also suitable for rotating instruments, such as used in podiatry

EAN: 4002516855194 / Material number: 12686260 / Old Material Number: 68549103D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
PWD 8534	•
Product category	
Insert/module for dental instruments	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of dental instruments	•
Capacity	
Sheaths [number]	48
Dimensions and weight	
External dimensions, gross height in mm	100
External dimensions, gross width in mm	265
External dimensions, gross depth in mm	180
Net weight in kg	0.30
Gross weight in kg	1.00