

PG 8597

Twin-tank water softener

for continuous provision of soft water with water hardness levels of up to 40 °dH

EAN: 4002516860372 / Material number: 12687200 / Old Material Number: 69859702D



- Output: continuous 19 l/min, max. 30 l/min
- Freestanding appliance with rollers, top loaded
- 1 container for 20 kg salt, 2 containers for 4.5 l resin
- Volume controlled two-chamber system, current-less operation
- The appliance is equipped with a water connection

Product category	
Tank dishwashers	•
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 8066	•
PG 8067	•
PG 8164	•
PG 8165	•
PG 8166	•
PG 8169	•
PG 8172	•
PG 8172 Eco	•
PTD 901	•
PTD 701	•
PTD 702	•
PTD 703	•
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PTD 704	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
Product category	
Water treatment	•
Product properties	
Material	Plastic
Colour	White
Dimensions and weight	
External dimensions, gross height in mm	680
External dimensions, gross width in mm	390
External dimensions, gross depth in mm	380
Net weight in kg	18.60
Gross weight in kg	19.00