

APWD 332
 Mesh tray
 flat, with 5 silicone dividers.



- Large perforation and almost total coverage during rinsing
- Ideal arrangement & securing of instruments in tray
- Convenient opening/closing with one hand
- Fits sterile supply container APWD 334
- Dimensions: 273 x 176 x 30 mm

EAN: 4002516910022 / Material number: 12752360 / Old Material Number: 69633200D

Product category	
Thermal disinfectors	•
Large-chamber thermal disinfectors, dental	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Small sterilisers	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PST 1710	•
PST 1720	•
PST 2210	•
PST 2220	•
PWD 7122	•
PWD 8682	•
PWD 8682 C	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Product properties	
Material	Stainless steel, Silicone
Colour	Stainless Steel, Yellow
Area of application	
Reprocessing of ophthalmology instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of VET instruments	•
Reprocessing of podiatry instruments	•
Reprocessing of surgical instruments	•
Reprocessing of lumened instruments	•
Dimensions and weight	
External dimensions, net height in mm	30
External dimensions, net width in mm	176
External dimensions, net depth in mm	273