

APWD 335 Sterile container for loading an APWD 333 mesh tray.



- Lid seal impermeable to germs
- Quick & contamination-free lid removal during sterile supply
- Rapid release filter for easy permanent filter changes
- Can be stacked – less space required
- Dimensions: 312 x 190 x 65 mm

EAN: 4002516909668 / Material number: 12752580 /
Old Material Number: 69633500D

Product category

Thermal disinfectors	•
Washer-disinfectors, medical	•
Small sterilisers	•

Machine affiliation

PST 1710	•
PST 1720	•
PST 2210	•
PST 2220	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•

Product category

Storage/transport system	•
--------------------------	---

Product properties

Material	Aluminium
Colour	Stainless Steel

Area of application

Reprocessing of ophthalmology instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of VET instruments	•
Reprocessing of podiatry instruments	•
Reprocessing of surgical instruments	•
Reprocessing of lumened instruments	•

Capacity

DIN mesh trays [number]	1
-------------------------	---

Dimensions and weight

External dimensions, net height in mm	65
External dimensions, net width in mm	190
External dimensions, net depth in mm	312