

APWD 336
Permanent filter
for up to 1000 cycles.



- Achieve good drying results in containers
- Resilient and durable
- 2 units per pack
- For use in APWD 334 and APWD 335

EAN: 4002516909903 / Material number: 12752660 /
Old Material Number: 69633600D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Small sterilisers	•
Machine affiliation	
PST 1710	•
PST 1720	•
PST 2210	•
PST 2220	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Product properties	
Material	Plastic
Colour	White
Area of application	
Reprocessing of ophthalmology instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of VET instruments	•
Reprocessing of podiatry instruments	•
Reprocessing of surgical instruments	•
Reprocessing of lumened instruments	•
Dimensions and weight	
External dimensions, net height in mm	215
External dimensions, net width in mm	94