

APWD 337

Holder

for identification plates in the size 56 x 18 mm.



EAN: 4002516910015 / Material number: 12752730 / Old Material Number: 69633700D

Product category	
Thermal disinfectors	•
Large-chamber thermal disinfectors, dental	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Small sterilisers	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PST 1710	•
PST 1720	•
PST 2210	•
PST 2220	•
PWD 7122	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of ophthalmology instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of VET instruments	•
Reprocessing of podiatry instruments	•
Reprocessing of surgical instruments	•
Reprocessing of lumened instruments	•
Dimensions and weight	
External dimensions, net height in mm	56
External dimensions, net width in mm	18