

ProCare Tex 11 COL - 500 g

Detergent for coloured items, powder, mildly alkaline, 500 g
for colour-fast cleaning of coloured items.



- For use in Miele Professional washing machines
- Highly concentrated detergent for coloured items
- Powerful grease and dirt removal that is gentle on materials and colours
- Full washing power even at 30 °C
- Very good protection against dye transfer
- Dermatologically tested

EAN: 4002516900085 / Material number: 12833710 / Old Material Number: 59710401V1

Product category	
Washing machines	•
Wash performance	
Removal of dirt	•
Removal of protein	•
Removal of grease	•
Antistatic action	•
Dispensing	
pH value at 20 °C	10,5
Manual dispensing	•
Can be used with all degrees of water hardness	•
How to use	
Self-service systems	•
Hotels and guesthouses	•
Spa, wellness and sports	•
Schools and nurseries	•
Small businesses	•
Retirement and care homes	•
Restaurants, bars and catering	•
Product properties	
Minimum storage temperature in °C	-20
Maximum storage temperature in °C	40
Dimensions and weight	
External dimensions, net height in mm	124
External dimensions, net width in mm	118
External dimensions, net depth in mm	118
External dimensions, gross height in mm	12.00
External dimensions, gross width in mm	300.00
External dimensions, gross depth in mm	400.00
Weight in g	500
Gross weight in kg	0.56
Product properties	
Medium type	Detergent for delicates/coloured items
Alkalinity	Mildly alkaline
Aggregate form	Powder
Composition	Surfactants, Fragrances, Enzymes, Sequestering agent
Biodegradable	•
Biologically degradable thanks to special-purpose enzymes	•
Phosphate-free	•
Particularly effective at low temperatures	•
Highly effective, even at 40 °C	•
Short activation time	•
High-level detergency combined with sparing use	•
Degrades without any residues	•
Free from micro-plastics	•
Simple handling	•

ProCare Tex 11 COL - 500 g

Detergent for coloured items, powder, mildly alkaline, 500 g
for colour-fast cleaning of coloured items.



EAN: 4002516900085 / Material number: 12833710 / Old Material Number: 59710401V1

Ideal additive for all detergents

•

Storage information

Minimum storage temperature in °C

-20

Maximum storage temperature in °C

40

Store in a dry location that is free of frost

•

Containers should always be re-sealed immediately once the required amount of product has been removed

•

Can be used after defrosting without loss of quality

•

Minimum shelf life in months

24

Dimensions and weight

External dimensions, net height in mm

124

External dimensions, net width in mm

118

External dimensions, net depth in mm

118

External dimensions, gross height in mm

12

External dimensions, gross width in mm

300

External dimensions, gross depth in mm

400

Weight in g

500

Gross weight in kg

0.56 kg

Compliance, test certificates, and symbols

GHS/CLP compliant

•