

ProCare Tex 70 S - 1 l

Fabric conditioner, liquid concentrate, acidic, 1 l
for care of fibres for long-term fabric touch and feel.



- For use in Miele Professional washing machines
- Ensures a particularly soft finish
- Good antistatic properties in the case of mixed fibre
- Excellent adsorption
- No encapsulated fragrances

EAN: 4002516900092 / Material number: 12833720 /
Old Material Number: 59704401V1

Product category	
Washing machines	•
Wash performance	
Antistatic action	•
Excellent adsorption	•
Dispensing	
Adjustment of dispensing concentration possible	•
pH value at 20 °C	4,5
Manual dispensing	•
Can be used with all degrees of water hardness	•
How to use	
Self-service systems	•
Hotels and guesthouses	•
Spa, wellness and sports	•
Schools and nurseries	•
Small businesses	•
Retirement and care homes	•
Restaurants, bars and catering	•
Product properties	
Minimum storage temperature in °C	4
Maximum storage temperature in °C	40
Dimensions and weight	
External dimensions, net height in mm	265
External dimensions, net width in mm	80
External dimensions, net depth in mm	80
External dimensions, gross height in mm	80.00
External dimensions, gross width in mm	80.00
External dimensions, gross depth in mm	270.00
Gross weight in kg	1.07
Fill volume in ml	1000
Product properties	
Medium type	Fabric conditioner
Alkalinity	Acidic
Aggregate form	Liquid concentrate
Composition	Fragrances
Phosphate-free	•
Free from micro-plastics	•
Ideal additive for all detergents	•
Storage information	
Minimum storage temperature in °C	4
Maximum storage temperature in °C	40
Store in a dry location that is free of frost	•
Containers should always be re-sealed immediately once the required amount of product has been removed	•
Dimensions and weight	
External dimensions, net height in mm	265
External dimensions, net width in mm	80
External dimensions, net depth in mm	80

ProCare Tex 70 S - 1 l

Fabric conditioner, liquid concentrate, acidic, 1 l
for care of fibres for long-term fabric touch and feel.



EAN: 4002516900092 / Material number: 12833720 /
Old Material Number: 59704401V1

External dimensions, gross height in mm	80
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	270
Gross weight in kg	1.07 kg
Fill volume in ml	1000

Compliance, test certificates, and symbols

GHS/CLP compliant

•