

Boost CX1 Blue Pulse
Bagless cylinder vacuum cleaners
for maximum extraction power in a compact design.



- Height-adjustable EasySlide telescopic tube
- Always stays on track thanks to TrackDrive
- Powerful cleaning performance thanks to vortex technology
- Versatile applications – universal floorhead
- Retains even the smallest particles – Hygiene AirClean filter

EAN: 4002516900481 / Material number: 12839360 /
Old Material Number: 41NRF340GB

Construction type	
Bagless vacuum cleaner	•
Product brand	Boost CX1
Motor type	PowerLine
Product type	SNRF3
Construction type and design	
Appliance colour	Obsidian black
Cleaning	
Vortex Technology	•
User convenience	
Electronic suction power control	4-stage rotary selector
Ergonomic carrying handle	•
Ergonomic handle	•
Park system	on one side
Mobility	
TrackDrive	•
Operating radius in m	10
Cable length in m	6.5
Filter system	
Exhaust filter	Hygiene AirClean filter
Technical data	
Power	890
Volume in dB	78.0
Width in mm	400
Height in mm	280
Depth in mm	28.00
Weight in kg	7.8
Nominal rating in W	890
Accessories included	
Upholstery nozzle	•
Crevice nozzle	•
Dusting brush	•
Floorhead SBD 355-3	•
Cleaning tool	No