

E 900-6 Mobile unit for the optimum cleaning of instruments on 6 levels.



- Equipped with connections for hollow instruments
- Loading capacity: 18 DIN mesh trays*
- Equipped with 3 integrated spray arms
- The loading areas for the individual levels are different
- Magnetic strip for mobile unit recognition (incl. magnets ML)

* DIN mesh tray 485 x 254 x 50 mm

EAN: 4002516932345 / Material number: 12889770 /
Old Material Number: 69590009D

Product category

- Large-capacity washer-disinfectors, laboratory
- Large-capacity washer-disinfectors, medical

Machine affiliation

- PG 8527
- PG 8528

Product category

- Mobile unit for inserts
- Mobile unit for mesh trays
- Mobile unit for MIS instruments

Area of application

- Reprocessing of laboratory glassware
- Reprocessing of surgical instruments

Product properties

- Material: Stainless steel
- Colour: Stainless Steel

Capacity

- DIN mesh trays [number]: 18
- SPRI mesh trays [number]: 12
- Lumened instruments [number]: 16

Dimensions and weight

- External dimensions, net height in mm: 553
- External dimensions, net width in mm: 636
- External dimensions, net depth in mm: 790
- External dimensions, gross height in mm: 680
- External dimensions, gross width in mm: 560
- External dimensions, gross depth in mm: 650
- Net weight in kg: 12.85
- Gross weight in kg: 12.85

Companion pack

- Magnetic strip for automatic mobile unit recognition (incl. magnets) [number]: 1