

E 913/3

Mobile unit

for the optimum loading of containers, lids and mesh trays on 3 levels.



- The loading areas for the individual levels are different
- Equipped with 2 integrated spray arms
- Equipped with connections for 5 hollow instruments
- Equipped with connection for hot air drying
- Magnetic strip for mobile unit recognition (incl. magnets ML)

EAN: 4002516949916 / Material number: 12890940 / Old Material Number: 69591306D

## Product category

Large-capacity washer-disinfectors, medical

•

## Machine affiliation

PG 8527

•

PG 8528

•

## Product category

Mobile unit for inserts

•

Mobile unit for containers

•

## Product properties

Material

Stainless steel

Colour

Stainless Steel

## Area of application

Reprocessing of surgical instruments

•

Reprocessing of containers

•

## Capacity

DIN mesh trays [number]

6

SPRI mesh trays [number]

2

DIN containers, height 15 cm [number]

3

Lumened instruments [number]

5

## Companion pack

Magnetic strip for automatic mobile unit recognition (incl. magnets) [number]

1