


E 940

Mobile unit

for the optimum loading of narrow necked glassware on 2 levels.



- The mobile injector unit can also be used without the upper level
- Magnetic strip for automatic mobile unit recognition
- Lower level equipped with a total of 115 nozzles with clips

EAN: 4002516948483 / Material number: 12895220 /
Old Material Number: 69594002D

Product category

Large-capacity washer-disinfectors, laboratory

•

Machine affiliation

PG 8527

•

PLW 8615

•

PLW 8616

•

PLW 8617

•

Product category

Mobile unit for narrow-necked glassware

•

Area of application

Reprocessing of laboratory glassware

•

Product properties

Material

Stainless steel, Plastic

Colour

Stainless Steel, Black

Capacity

Narrow-necked glassware [number]

115

Companion pack

Magnetic strip for automatic mobile unit recognition (incl. magnets) [number]

1