

E 941
Mobile unit
for modules on 2 levels.



- Suitable for up to 2 modules per level depending on the size of the load
- Coupling of the mobile unit to the water and drying air supply
- The clearances for both levels are different

EAN: 4002516949404 / Material number: 12895230 /
Old Material Number: 69594102D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Mobile unit for orthopaedic instruments	•
Other mobile units	•
Area of application	
Reprocessing of laboratory glassware	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Companion pack	
Magnetic strip for automatic mobile unit recognition (incl. magnets) [number]	1