

ProCare Universal 70 ST - 750 ml
Multi-purpose cleaner, mildly alkaline, 750 ml
for gently removing grease residues and contamination.



- Effective high concentrate for optimum cleaning results
- Excellent dissolving power for oil and grease
- Very gentle on materials
- Recommended for Miele Professional and Steelco stainless steel machines
- Awarded EU Ecolabel
- Recyclable packaging

EAN: 4002516954507 / Material number: 12938290 /
Old Material Number: 66917702V1

Product category	
Fresh water dishwashers	•
ProfiLine dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Washing machines	•
Tumble dryers	•
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Small sterilisers	•
Applications	
Cleaning of stainless steel surfaces	•
Manual care and maintenance	•
Removal of starch deposits	•
Removal of stubborn or dried-on residues	•
Removal of grease	•
Removal of oil	•
Removal of protein	•
Removal of blood	•
Dispensing	
pH value at 20 °C	10,8
Product properties	
Medium type	Stainless steel cleaner
Aggregate form	Liquid concentrate
Ingredients	Glycol ethers, Sequestering agent (no phosphates), Solubility promoter, Anionic surfactants, Non-ionic surfactants
Minimum storage temperature in °C	5
Maximum storage temperature in °C	25
Recyclable packaging	•
Good cleaning capability	•
Does not contain phosphates	•
Does not contain phosphorus compounds	•
Does not contain silicates	•
Very gentle on materials	•
Intensive cleaning agent containing surfactants	•
Strong emulsifying and dispersing effect	•
Very good material compatibility	•
Dimensions and weight	
External dimensions, net height in mm	297
External dimensions, net width in mm	94
External dimensions, net depth in mm	63
External dimensions, gross height in mm	60.00
External dimensions, gross width in mm	90.00
External dimensions, gross depth in mm	260.00

ProCare Universal 70 ST - 750 ml

Multi-purpose cleaner, mildly alkaline, 750 ml

for gently removing grease residues and contamination.



EAN: 4002516954507 / Material number: 12938290 /

Old Material Number: 66917702V1

Gross weight in kg	0.81
--------------------	------

Fill volume in ml	750
-------------------	-----

Compliance, test certificates, and symbols	
---	--

EU eco label	•
--------------	---

Verification numbers	
-----------------------------	--

EU Eco label: DE/020/507	•
--------------------------	---