

E 315/1

Full demineralisation single use resin

For the regeneration of full demineralisation cartridges

EAN: 4002516957515 / Material number: 12945010 / Old Material Number: 69431502D



- For use in the LP 2800/1 empty cartridge
- 20 l of homogeneously mixed resins
- Carton with 2 x 10 l bags
- Vacuum sealed in plastic bags
- Including replacement filter bag

Product category	
Fresh-water dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Laboratory washers	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8067	•
PG 8041	•
PG 8043	•
PG 8045	•
PG 8055	•
PG 8055 U	•
PG 8056	•
PG 8056 U	•
PG 8057 TD	•
PG 8057 TD U	•
PG 8058	•
PG 8058 U	•
PG 8059	•
PG 8059 U	•
PG 8060	•
PG 8061	•
PG 8061 U	•
PG 8063	•
PG 8164	•
PG 8165	•
PG 8166	•
PG 8169	•
PG 8172	•
PG 8172 Eco	•
PTD 901	•
PG 8096	•
PG 8096 U	•
PG 8099	•
PG 8099 U	•
PFD 400	•
PFD 400 U	•
PFD 401	•
PFD 401 U	•
PFD 402	•
PFD 402 U	•
PFD 404	•
PFD 404 U	•
PFD 405	•
PFD 405 U	•
PFD 407	•

E 315/1

Full demineralisation single use resin

For the regeneration of full demineralisation cartridges

EAN: 4002516957515 / Material number: 12945010 / Old Material Number: 69431502D



PFD 407 U	•
PFD 408 U	•
PTD 701	•
PTD 702	•
PTD 703	•
PTD 704	•
PLW 8683	•
PLW 8683 CD	•
PWD 8682	•
PWD 8682 C	•
PWD 8692	•
Product category	
Water treatment	•
Dimensions and weight	
External dimensions, gross height in mm	285
External dimensions, gross width in mm	260
External dimensions, gross depth in mm	440
Net weight in kg	13.77
Gross weight in kg	14.27