

ProCare Lab 10 AP - 5 l

liquid cleaning agent, alkaline, 5 l

for machine-based universal cleaning of laboratory glassware and utensils.



- Alkaline cleaning agent with a high percentage of dispersants
- Free from surfactants, phosphates and oxidants
- Specially developed for Miele Professional washer-disinfectors
- Suitable for stainless steel, laboratory glassware, ceramics and plastics
- For cleaning stubborn or dried-on residues

EAN: 4029577160142 / Material number: 12958900 / Old Material Number: 62924102V1

Product category	
Laboratory washers	•
Large-chamber laboratory washers	•
Applications	
Removal of stubborn or dried-on residues	•
Removal of grease	•
Removal of oil	•
Removal of organic residues	•
Dispensing	
Dispensing recommendation per wash cycle	2-4 ml/l (0.2-0.4%)
Large opening for suction lance	•
Automatic dispensing	•
Recommended water hardness level in mmol/l	0-0.5
Main applications	
Medical laboratories	•
Biological laboratories	•
Microbiological laboratories	•
Chemical laboratories	•
Water laboratories	•
Laboratories for the food industry	•
Laboratories for the cosmetics industry	•
Laboratories for the mineral oil industry	•
Industrial laboratories	•
Phosphate industry	•
Colour industry	•
Laboratories for the pharmaceutical industry	•
Product properties	
Medium type	Alkaline cleaning agent
Aggregate form	Liquid concentrate
Ingredients	Potassium hydroxide, Polycarboxylates
Minimum storage temperature in °C	0
Maximum storage temperature in °C	30
pH value at 20 °C	14,0
Suitable for plastic items	•
Suitable for stainless steel	•
Suitable for laboratory glassware	•
Suitable for ceramic items	•
Not suitable for aluminium	•
Not suitable for light and non-ferrous metals	•
Does not contain oxidants	•
Does not contain phosphates	•
Does not contain surfactants	•
Free from microplastics	•
Good cleaning capability	•
High percentage of dispersants	•
Dimensions and weight	
External dimensions, net height in mm	217

ProCare Lab 10 AP - 5 l
liquid cleaning agent, alkaline, 5 l
for machine-based universal cleaning of laboratory glassware and utensils.



EAN: 4029577160142 / Material number: 12958900 / Old Material Number: 62924102V1

External dimensions, net width in mm	142
External dimensions, net depth in mm	248
External dimensions, gross height in mm	250
External dimensions, gross width in mm	190
External dimensions, gross depth in mm	150
Gross weight in kg	7.50
Fill volume in ml	5000
Compliance, test certificates, and symbols	
GHS/CLP compliant	•