

ProCare Lab 10 MA - 5 l

liquid cleaning agent, alkaline, 5 l

for material-saving machine cleaning of laboratory glassware and utensils.



- Universal alkaline cleaning agent
- Free from oxidants and surfactants
- Specially developed for Miele Professional washer-disinfectors
- Suitable for stainless steel, laboratory glassware, ceramics and plastics
- Effective cleaning for light to medium soiling

EAN: 4029577160180 / Material number: 12959990 / Old Material Number: 62923102V1

Product category	
Laboratory washers	•
Large-chamber laboratory washers	•
Applications	
Removal of protein	•
Removal of blood	•
Removal of residues of nutrient mediums	•
Removal of organic residues	•
Dispensing	
Dispensing recommendation per wash cycle	2-4 ml/l (0.2-0.4%)
Large opening for suction lance	•
Automatic dispensing	•
Suitable for all levels of water hardness	•
Main applications	
Medical laboratories	•
Biological laboratories	•
Chemical laboratories	•
Water laboratories	•
Laboratories for the food industry	•
Laboratories for the cosmetics industry	•
Serological institutes	•
Laboratories for the pharmaceutical industry	•
Product properties	
Medium type	Alkaline cleaning agent
Aggregate form	Liquid concentrate
Ingredients	Potassium hydroxide, Phosphate, Silicates
Minimum storage temperature in °C	-15
Maximum storage temperature in °C	30
pH value at 20 °C	14,0
Suitable for plastic items	•
Suitable for stainless steel	•
Suitable for laboratory glassware	•
Suitable for ceramic items	•
Limited suitability for non-ferrous metals	•
Limited suitability for aluminium	•
Not suitable for aluminium	•
Not suitable for light alloys	•
Not suitable for light metals	•
Does not contain oxidants	•
Does not contain surfactants	•
Free from microplastics	•
Good cleaning capability	•
Very gentle on materials	•
Wide range of uses	•
Dimensions and weight	
External dimensions, net height in mm	217
External dimensions, net width in mm	142

ProCare Lab 10 MA - 5 l
liquid cleaning agent, alkaline, 5 l
for material-saving machine cleaning of laboratory glassware and utensils.



EAN: 4029577160180 / Material number: 12959990 / Old Material Number: 62923102V1

External dimensions, net depth in mm	248
External dimensions, gross height in mm	250
External dimensions, gross width in mm	190
External dimensions, gross depth in mm	150
Gross weight in kg	8.00
Fill volume in ml	5000
Compliance, test certificates, and symbols	
GHS/CLP compliant	•