

ProCare Lab 10 MA - 10 l

liquid cleaning agent, alkaline, 10 l

for material-saving machine cleaning of laboratory glassware and utensils.

EAN: 4029577160197 / Material number: 12960010 / Old Material Number: 62923202V1



- Universal alkaline cleaning agent
- Free from oxidants and surfactants
- Specially developed for Miele Professional washer-disinfectors
- Suitable for stainless steel, laboratory glassware, ceramics and plastics
- Effective cleaning for light to medium soiling

Product category	
Laboratory washers	•
Large-chamber laboratory washers	•
Applications	
Removal of protein	•
Removal of blood	•
Removal of residues of nutrient mediums	•
Removal of organic residues	•
Dispensing	
Dispensing recommendation per wash cycle	2-4 ml/l (0.2-0.4%)
Large opening for suction lance	•
Automatic dispensing	•
Suitable for all levels of water hardness	•
Main applications	
Medical laboratories	•
Biological laboratories	•
Chemical laboratories	•
Water laboratories	•
Laboratories for the food industry	•
Laboratories for the cosmetics industry	•
Serological institutes	•
Laboratories for the pharmaceutical industry	•
Product properties	
Medium type	Alkaline cleaning agent
Aggregate form	Liquid concentrate
Ingredients	Potassium hydroxide, Phosphate, Silicates
Minimum storage temperature in °C	-15
Maximum storage temperature in °C	30
pH value at 20 °C	14,0
Suitable for plastic items	•
Suitable for stainless steel	•
Suitable for laboratory glassware	•
Suitable for ceramic items	•
Limited suitability for non-ferrous metals	•
Limited suitability for aluminium	•
Not suitable for aluminium	•
Not suitable for light alloys	•
Not suitable for light metals	•
Does not contain oxidants	•
Does not contain surfactants	•
Free from microplastics	•
Good cleaning capability	•
Very gentle on materials	•
Wide range of uses	•
Dimensions and weight	
External dimensions, net height in mm	308
External dimensions, net width in mm	191

ProCare Lab 10 MA - 10 l

liquid cleaning agent, alkaline, 10 l

for material-saving machine cleaning of laboratory glassware and utensils.

EAN: 4029577160197 / Material number: 12960010 / Old Material Number: 62923202V1



External dimensions, net depth in mm	234
External dimensions, gross height in mm	320
External dimensions, gross width in mm	230
External dimensions, gross depth in mm	190
Gross weight in kg	15.50
Fill volume in ml	10000

Compliance, test certificates, and symbols

GHS/CLP compliant

•