

ProCare Lab 30 C - 5 l

Neutralising agent based on citric acid, acidic, 5 l

for material-saving machine cleaning of laboratory glassware and utensils.



- To be used as neutralising agent or acid pre-cleaner
- Free from surfactants and phosphates
- Specially developed for Miele Professional washer-disinfectors
- Suitable for stainless steel, laboratory glassware, ceramics and plastics
- Removes residues that are soluble in acids from sensitive materials

EAN: 4029577160227 / Material number: 12960630 / Old Material Number: 62927102V1

Product category	
Laboratory washers	•
Large-chamber laboratory washers	•
Applications	
Removal of limescale and other acid-soluble impurities	•
Dispensing	
Dispensing recommendation per wash cycle	1–5 ml/l (0.1–0.5 %)
Large opening for suction lance	•
Automatic dispensing	•
Suitable for all levels of water hardness	•
Main applications	
Medical laboratories	•
Biological laboratories	•
Microbiological laboratories	•
Chemical laboratories	•
Water laboratories	•
Laboratories for the food industry	•
Laboratories for the cosmetics industry	•
Laboratories for the phosphate industry	•
Laboratories for the pharmaceutical industry	•
Product properties	
Medium type	Acidic neutralising agent
Aggregate form	Liquid concentrate
Ingredients	Citric acid, Preservatives
Minimum storage temperature in °C	-3
Maximum storage temperature in °C	30
pH value at 20 °C	1,0
Suitable for plastic items	•
Suitable for stainless steel	•
Suitable for laboratory glassware	•
Suitable for ceramic items	•
Suitable for sensitive materials	•
Suitable for neutralisation of non-anodised aluminium	•
Not suitable for anodised aluminium	•
Not suitable for light alloys	•
Not suitable for acidic pre-cleaning of non-anodised aluminium	•
Does not contain phosphates	•
Does not contain surfactants	•
Free from microplastics	•
Very gentle on materials	•
Acidic pre-cleaning	•
Neutralisation of carried-over alkali residues	•
Dimensions and weight	
External dimensions, net height in mm	217
External dimensions, net width in mm	142
External dimensions, net depth in mm	248
External dimensions, gross height in mm	250

ProCare Lab 30 C - 5 l
Neutralising agent based on citric acid, acidic, 5 l
for material-saving machine cleaning of laboratory glassware and utensils.



EAN: 4029577160227 / Material number: 12960630 / Old Material Number: 62927102V1

External dimensions, gross width in mm	190
External dimensions, gross depth in mm	150
Gross weight in kg	6.50
Fill volume in ml	5000
Compliance, test certificates, and symbols	
GHS/CLP compliant	•