

ProCare Lab 30 P - 5 l

Neutralising agent based on phosphoric acid, acidic, 5 l
for machine cleaning of laboratory glassware and utensils.



- To be used as neutralising agent or acid pre-cleaner
- Does not contain surfactants
- Specially developed for Miele Professional washer-disinfectors
- Suitable for stainless steel, laboratory glassware, ceramics and plastics
- Removes stubborn residues that are soluble in acids

EAN: 4029577160241 / Material number: 12960700 / Old Material Number: 62928102V1

Product category	
Laboratory washers	•
Large-chamber laboratory washers	•
Applications	
Removal of limescale and other acid-soluble impurities	•
Dispensing	
Dispensing recommendation per wash cycle	1-4 ml/l (0.1-0.4%)
Large opening for suction lance	•
Automatic dispensing	•
Suitable for all levels of water hardness	•
Main applications	
Medical laboratories	•
Biological laboratories	•
Microbiological laboratories	•
Chemical laboratories	•
Laboratories for the food industry	•
Laboratories for the cosmetics industry	•
Laboratories for the mineral oil industry	•
Industrial laboratories	•
Laboratories for the pharmaceutical industry	•
Product properties	
Medium type	Acidic neutralising agent
Aggregate form	Liquid concentrate
Ingredients	Phosphoric acid, Citric acid
Minimum storage temperature in °C	-15
Maximum storage temperature in °C	30
pH value at 20 °C	0,7
Suitable for plastic items	•
Suitable for stainless steel	•
Suitable for laboratory glassware	•
Suitable for ceramic items	•
Not suitable for aluminium	•
Not suitable for light and non-ferrous metals	•
Not suitable for light alloys	•
Does not contain surfactants	•
Free from microplastics	•
Acidic pre-cleaning	•
Excellent rinsing properties	•
Neutralisation of carried-over alkali residues	•
Dimensions and weight	
External dimensions, net height in mm	217
External dimensions, net width in mm	142
External dimensions, net depth in mm	248
External dimensions, gross height in mm	200
External dimensions, gross width in mm	140
External dimensions, gross depth in mm	240
Gross weight in kg	8.00
Fill volume in ml	5000

ProCare Lab 30 P - 5 l

Neutralising agent based on phosphoric acid, acidic, 5 l
for machine cleaning of laboratory glassware and utensils.



EAN: 4029577160241 / Material number: 12960700 / Old Material Number: 62928102V1

Compliance, test certificates, and symbols

GHS/CLP compliant

•