

ProCare Universal 77 A - 500 g
Machine cleaner, alkaline, 500 g
for removing persistent starch deposits.



- Removal of grease deposits
- Reduction in odours
- Specially developed for Miele Professional machines

EAN: 4002517008070 / Material number: 13068710 / Old Material Number: 66921702V1

Product category	
Fresh water dishwashers	•
ProfiLine dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Washing machines	•
Laboratory washers	•
Large-chamber laboratory washers	•
Thermal disinfectors	•
Large-chamber thermal disinfectors	•
Washer-disinfectors, dental	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Applications	
Removal of starch deposits	•
Removal of grease	•
Removal of oil	•
Removal of protein	•
Dispensing	
Dispensing recommendation per wash cycle	0.5–1 dose
Suitable for all levels of water hardness	•
Manual dispensing	•
Product properties	
Medium type	Machine cleaner
Aggregate form	Powder
Ingredients	Oxygen-based bleaching agent, Polycarboxylates, Surfactants
Minimum storage temperature in °C	5
Maximum storage temperature in °C	35
pH value at 20 °C	11,0
Good cleaning capability	•
Does not contain phosphates	•
Very gentle on materials	•
High percentage of dispersants	•
Intensive cleaning agent containing surfactants	•
Free from fragrances	•
Free from microplastics	•
Dimensions and weight	
External dimensions, net height in mm	133
External dimensions, net width in mm	80
External dimensions, net depth in mm	80
External dimensions, gross height in mm	130
External dimensions, gross width in mm	80
External dimensions, gross depth in mm	80
Weight in g	500
Gross weight in kg	0.57
Compliance, test certificates, and symbols	
GHS/CLP compliant	•